



Recombinant *Synechococcus* sp. Uncharacterized protein SYNPC7002_A1392 (SYNPC7002_A1392)

Product Code	CSB-MP329128FQA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P31972
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus</i> sp. (strain ATCC 27264 / PCC 7002 / PR-6) (<i>Agmenellum quadruplicatum</i>)
Purity	>85% (SDS-PAGE)
Sequence	MTAHNFFWDG DWFSPQASQD QNLIAQLDFV PGLRELLFIR QAHALEHGTV WVLTELAERY PTAWRHHYAE LSGMSTGQGF YLFGGVDKTD LYRAVSTAGD RLRSGEWNLA IHPRCGTNIS VTLLLTAGLA GGAAWLLPKD PLTQLFGMGT AAATAWAIAP DLGKYAQQHI TTAIPLNLAL EEIQENTDES GNPSHFVRLR WQDAQ
Source	Mammalian cell
Target Names	SYNPC7002_A1392
Protein Names	Recommended name: Uncharacterized protein SYNPC7002_A1392
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.