



Recombinant Human Protein FAM181B (FAM181B)

Product Code	CSB-MP008152HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NEQ2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAVQAALLST HPFVPGFGG SPDGLGGAFG ALDKGCCFED DETGAPAGAL LSGAEGGDVR EATRDLLSFI DSASSNIKLA LDKPGKSKRK VNHRKYLQKQ IKRCSGLMGA APPGPPPSA ADTPAKRPLA APSAPTVAAP AHGKAAPRRE ASQAAAAASL QSRSLAALFD SLRHVPGGAE PAGGEVAAPA AGLGGAGTGG AGGDVAGPAG ATAIPGARKV PLRARNLPPS FFTEPSRAGG GGCGPSGPDV SLGDLEKGAE AVEFFELLGP DYGAGTEAAV LLAAEPLDVF PAGASVLRGP PELEPGLFEP PPAVVGNNLY PEPWSVPGCS PTKKSPLTAP RGGLTLNEPL SPLYPAAADS PGGEDGRGHL ASFAPFFPDC ALPPPPPHQ VSYDYSAGYS RTAYSSLWRS DGVWEGAPGE EGAHRD
Source	Mammalian cell
Target Names	FAM181B
Protein Names	Recommended name: Protein FAM181B
Expression Region	1-426
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.