



Recombinant Escherichia coli Uncharacterized protein yjcF (yjcF)

Product Code	CSB-YP329193ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32704
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MRYNGLNNMF FPLCLINDNH SVTSPSHTKK TKSDNYSKHH KNTLIDNKAL SLFKMDDHEK VIGLIQKMKR IYDSLPSGKI TKETDRKIHK YFIDIASHAN NKDDRITRR VYLNKDKEVS IKVVYFINNV TVHNNTIEIP QTVNGGYDFS HLSLKGIVIK DEDLSNSNFA GCRLQNAIFQ DCNMYKTNFN FAIMEKILFD NCILDDSNA QIKMTDGTLN SCSAMHVQFY NATMNRANIK NTFLDYSNFY MAYMAEVNLY KVIAPYINLF RADLSFSKLD LINFEHADLS RVNLNKATLQ NINLIDSKLF FTSLTNTFLE MVICTDSNMA NVNFNNANLS NCHFNCVLT KAWMFNIRLY RVNFDEASVQ GMGITLRGE ENISINS DIL VTLQKFFFEED CATHTGMSQT EDNLHAVAMK ITADIMQDAD
Source	Yeast
Target Names	yjcF
Protein Names	Recommended name: Uncharacterized protein yjcF
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.