



Recombinant *Saccharomyces cerevisiae* S-(hydroxymethyl)glutathione dehydrogenase

Product Code	CSB-BP329196SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32771
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSAATVGKPI KCIAAVAYDA KKPLSVEEIT VDAKPAHEVR IKIEYTAVCH TDAYTLGSD PEGLFPCVLG HEGAGIVESV GDDVITVKPG DHVIALYTAE CGKCKFCTSG KTNLCGAVRA TQGKGVMPDG TTRFHNAKGE DIYHFMGCST FSEYTVVADV SVVAIDPKAP LDAACLLGCG VTTGFGAALK TANVQKGDV AVFGCGTVGL SVIQGAKLRG ASKIIAIDIN NKKKQYCSQF GATDFVNPKE DLAKDQTIVE KLIEMTDGGL DFTFDCTGNT KIMRDALEAC HKGWGQSIII GVAAAGEEIS TRPFQLVTGR VWKGSFAFGGI KGRSEMGGI KDYQKGALKV EEFITHRRPF KEINQAFEDL HNGDCLRTVL KSDEIK
Source	Baculovirus
Target Names	SFA1
Protein Names	Recommended name: S-(hydroxymethyl)glutathione dehydrogenase EC= 1.1.1.284 Alternative name(s): Alcohol dehydrogenase SFA EC= 1.1.1.1 Glutathione-dependent formaldehyde dehydrogenase Short name= FALDH Short name= FD
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.