



Recombinant *Saccharomyces cerevisiae* Methionyl-tRNA formyltransferase, mitochondrial (FMT1)

Product Code	CSB-MP329198SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32785
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	VQPL NVLFFGSDTF SNFSLQALNE LRQNGSCGI VDNIQVVTRS PKWCGRQKSI LKYPPIFDMA EKLQLRPIT CDTKQEMLAL SKLTPSRQGN PENDGSGAPF NAIIVSFGK LIPGDLIRAV PLALNVHPSL LPRHKGSAPI QRALLEGDTY TGVTIQTLHP DRFDHGAIVA QTEPLAIATM LSKGRVNDST ADFNSEGLPR RTAILMDQLG ALGAQLLGQT LRERLYLPQN RVQAPTAYKP SYAHRITTED KRIHWARDSA AELLNKLETL GPLHAFKEAT AARKDAQNSV LKRIKFHECK VMRDARLDNG SKPGMFKYDD IKDCILVTCR GNLLLCVSRL QFEGFAVERA GQFMARLRKR CGALSEKLVF L
Source	Mammalian cell
Target Names	FMT1
Protein Names	Recommended name: Methionyl-tRNA formyltransferase, mitochondrial Short name= MtFMT EC= 2.1.2.9
Expression Region	27-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.