



# Recombinant *Saccharomyces cerevisiae* Chromatin structure-remodeling complex subunit RSC7 (NPL6)

<b>Product Code</b>	CSB-EP329205SVG-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P32832
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSDSEGLAS EVEHEKRSRS TSNRPNYAID TEDLDIDEND ENEDDDYREE EANEGVNEEE ISDEEEQINK SGRNKRRHVD EEEDLSEDKG VTRSRNRSKF KKPVFPGIDD AEENLNPLKV VNEEYVLPDD PEGETKITAD GDLLGGREFL VRTFTLTEKG NRKFMLATEP ARIVGFRDSY LFFQTHPNLY KFILNQTQKN DLIDRGVLPY SYRNRQIALV TARGVFKEFG AKIIRGGKHI TDDYYASELR TKGNVIEGKL AGDPIDKSAR ALETMMYPAS ENGINPAKNQ VEFFEHRPHG HMSNSNIAS GSKLSSTNWL YQHSAACSRF NSDLFYDRVK VLLVDQQGLR DAYTNILHIP ESTQSTTVLG WRRSKNDSPS DTSIVYETVI HDNDLNKPKT GLSEIPKEIY EDVVDEDVLR AITEQQNF EK CNEYI
<b>Source</b>	E.coli
<b>Target Names</b>	NPL6
<b>Protein Names</b>	Recommended name: Chromatin structure-remodeling complex subunit RSC7 Alternative name(s): Nuclear protein localization protein 6 Remodel the structure of chromatin complex subunit 7
<b>Expression Region</b>	1-435
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.