



# Recombinant Escherichia coli Putative acyl-CoA dehydrogenase AidB (aidB)

<b>Product Code</b>	CSB-BP329241ENV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P33224
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MHWQHTTVFN QPIPLNNSNL YLSDGALCEA VTREGAGWDS DFLASIGQQL GTAESLELGR LANVNPELL RYDAQGRRLD DVRFHPAWHL LMQALCTNRV HNLAWEDAR SGAFVARAAR FMLHAQVEAG SLCPITMTFA ATPLLLQMLP APFQDWTTP LLSDRYDSHLL PGGQKRGLLI GMGMTEKQGG SDVMSNTTRA ERLEDGSYRL VGHKWFFSVP QSDAHLVLAQ TAGGLSCFFV PRFLPDGQRN AIRLERLKDK LGNRSNASCE VEFQDAIGWL LGLEGE GIRL ILKMGGMTRF DCALGSHAMM RRAFSLAIYH AHQRHVFVGNP LIQQPLMRHV LSRMALQLEG QTALLFRLAR AWDRRADAKE ALWARLFTPA AKFVICKRGM PFVAEAMEVL GGIGYCEESE LPRLYREMPV NSIWEGSGNI MCLDVLRLV LN KQAGVYD LLS EAFVEVKGQD RYFDRAVRR LQQQLRKPAEE LGREITHQLF LLGCGAQMLK YASPPMAQAW CQVMLDTRGG VRLSEIQIQND LLLRATGGVC V
<b>Source</b>	Baculovirus
<b>Target Names</b>	aidB
<b>Protein Names</b>	Recommended name: Putative acyl-CoA dehydrogenase AidB EC= 1.3.99.-
<b>Expression Region</b>	1-541
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.