



# Recombinant Moraxella bovis Type-2 restriction enzyme MbolI (mbolIR)

<b>Product Code</b>	CSB-MP329382MHZ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P23191
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Moraxella bovis
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKNYVSNINL GNSSLKFIDE RLQSENYRGI HLSQHNRYDL PKLIDILTLL NKHAPNQSLM QIRTTDISKR PQNIPEEQSY AEFCEAKSL TNIGTQDAMR KNLFVDFARM GLINRYNDKK VLTDPFKRGV TKYVALSDMG VKLIDPKLDI LSKNLIFSLS LNKLLTGFVE DVLSLLTNSD LKEISFDEFM LFVSAMNCNF NFSISTEQCE SLIKEYRLLS RVQKNAVIDT LKSELIPDNF NGDKKDKRDY HNWANENQQI WTLFENIPFF IMEKDSRCLI LITSDVDLSK YSKSKMKRSQ QAKNDYFKHH KVNKIKGYEL DHIIPLEAE SVDEYRYLDN WLNLLYIDGK THAIKSQSGS KYYIFTFDDN DYNQIYFLDT QGDKLSINND DTALFDKNKV PKIYEYNQNF INAKTS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	mbolIR
<b>Protein Names</b>	Recommended name: Type-2 restriction enzyme MbolI Short name= R.MbolI EC= 3.1.21.4Alternative name(s): Endonuclease MbolI Type II restriction enzyme MbolI Type IIS restriction enzyme MbolI
<b>Expression Region</b>	1-416
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.