



Recombinant Human WAS/WASL-interacting protein family member 3 (WIPF3)

Product Code	CSB-BP026116HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NGB9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPVPPPPPPP LPPPPPPLGA PPPPPPSAPP VSTDTSSLRR ADPKGRSALL ADIQQGTRLR KVTQINDRSA PQIESSKGTN KEGGGSANTR GASTPPTLGD LFAGGFPVLR PAGQRDVAGG KTGQGGPSRA PSPRLPNKTI SGPLIPPASP RLGNTSEAHG AARTAPPRPN VPAPPPPTPP PPPPPLPPPL PSSSPIKTPL VSPPGPLTKG NLPVVAPPVP CAPPPIPPPPP PPTPPPLPPA SVLSDKAVKP QLAPLHLPPI PPPLPLLPPC GYPGLKAEPA SPAQDAQEPP APPPPLPPYA SCSPRASLPA PPLPGVNSSS ETPPPLPPKS PSFQAPPQKA GAQALPAPPA PPGSQPFLQK KRHGRPGAGG GKLNPPPAPP ARSPTTELSS KSQQATAWTP TQQPGGQLRN GSLHIIDDFE SKFTFHSVED FPPPDEYKPC QKIYPSKIPR SRTPGPWLQA EAVGQSSDDI KGRNSQLSLK TLR
Source	Baculovirus
Target Names	WIPF3
Protein Names	Recommended name: WAS/WASL-interacting protein family member 3 Alternative name(s): Corticosteroids and regional expression protein 16 homolog
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.