



# Recombinant Human WAS/WASL-interacting protein family member 3 (WIPF3)

<b>Product Code</b>	CSB-EP026116HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A6NGB9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPVPPPPPPP LPPPPPPLGA PPPPPPSAPP VSTDTSLLRR ADPKGRSALL ADIQQGTRLR KVTQINDRSA PQIESSKGTN KEGGGSANTR GASTPPTLGD LFAGGFVLR PAGQRDVAGG KTGQGGPSRA PSPRLPNKTI SGPLIPPASP RLGNTSEAHG AARTAPPRPN VPAPPPPTPP PPPPPLPPPL PSSSPIKTPL VSPPGPLTKG NLPVVAPPVP CAPPPIPPPP PPTPPPLPPA SVLSDKAVKP QLAPLHLPI PPPLPLLPPC GYPGLKAEPA SPAQDAQEPP APPPPLPPYA SCSPRASLPA PPLPGVNSSS ETPPPLPPKS PSFQAPPQKA GAQALPAPPA PPGSQPFLQK KRHGRPGAGG GKLNPPAPP ARSPTELSS KSQQATAWTP TQQPGQLRN GSLHIIDDFE SKFTFHSVED FPPPDEYKPC QKIYPSKIPR SRTPGPWLQA EAVGQSSDDI KGRNSQLSLK TLR
<b>Source</b>	E.coli
<b>Target Names</b>	WIPF3
<b>Protein Names</b>	Recommended name: WAS/WASL-interacting protein family member 3 Alternative name(s): Corticosteroids and regional expression protein 16 homolog
<b>Expression Region</b>	1-483
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.