



Recombinant *Xenopus laevis* Guanine nucleotide-binding protein G (s) subunit alpha

Product Code	CSB-BP329505XBE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P24799
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	GCLGNSKTE DQRNEEKVQR ETNKKIEKQL QKDKQVYRAT HRLLLLGAGE SGKSSIVKQM RILHVNGFNA EEKTKVQDI KNNIKEAIET IVTAMGNLSP PVELVNPENQ FRIDYILNLP NYKDFEFSPE FYEHTKTLWQ DEGVRACYER SNEYQLIDCA QYFLDKIDIV KQNDYTPSDQ DLLRCRVLTS GIFETKFQVD KVNFMFDVG GQRDERRKWI QCFNDVTAIL FVASSSYNM VIREDNHTNR LQEALNLFKS IWNNRWLRTI SVILFLNKQD LLAEKVNAGK SKIEDYFPEF ARYTTPDDAT PEVGEDPRVT RAKYFIRDEF LRISTASGDG RHYCYPHFTC AVDTENIRRV FNDCRDIIQR MHLRQYELL
Source	Baculovirus
Target Names	gnas
Protein Names	Recommended name: Guanine nucleotide-binding protein G(s) subunit alpha Alternative name(s): Adenylate cyclase-stimulating G alpha protein
Expression Region	2-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.