



Recombinant *Saccharomyces cerevisiae* Putative uncharacterized protein YCR013C (YCR013C)

Product Code	CSB-EP329563SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P25614
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MGKEKRKKLI YRFQFNSIYF FSDKKATPGN SLPSNNSKEA PPPVETWEIL SVTPYFLATV AVSPPPMMTV LPAAELLTTS SNKALVPAAN FSNSKTPGGP FQTMVLAFAT VAANNFLDSG PLSNPCQPAG IPSLSVTVLV LASAEKASAM MKSTGKTTST PLALAFSINF GTISAPALSK MESPISVFSK TFLKVKAIPP PMMIESTLSN KLSIN
Source	E.coli
Target Names	YCR013C
Protein Names	Recommended name: Putative uncharacterized protein YCR013C
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.