



Recombinant Bovine coronavirus Nucleoprotein (N)

Product Code	CSB-YP329597BJL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P26020
Product Type	Recombinant Protein
Immunogen Species	Bovine coronavirus (strain OK-0514) (BCoV) (BCV)
Purity	>85% (SDS-PAGE)
Sequence	MSFTPGKQSS SRASSGNRSG NGILKWADQS DQSRNVQTRG RRAQPKQTAT SQQPSGGNVV PYYSWFSGIT QFQKGKEFEF AEGQGVPIAP GVPATEAKGY WYRHNRRSFK TADGNQRQLL PRWYFYLLGT GPHAKDQYGT DIDGVFVVAS NQADVNTPAD ILDRDPSSDE AIPTRFPPGT VLPQGYIIEG SGRSAPNSRS TSRASSRASS AGRSRANSR NRTPTSGVTP DMADQIASLV LAKLGKDATK PQQVTKQTAK EIRQKILNKP RQKRSPNKQC TVQQCFGKRG PNQNFGGGEM LKLGTSDPQF PILAELAPTA GAFFFGSRLE LAKVQNLSGN LDEPQKDVEE LRYNGAIRFD STLSGFETIM KVLNENLNAY QQQDGMMNMS PKPQRQRGQK NGQGENDNIS VAAPKSRVQQ NKSRELTAEI ISLLKKMDEP YTEDTSEI
Source	Yeast
Target Names	N
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein Short name= NC Short name= Protein N
Expression Region	1-448
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.