



Recombinant Tomato yellow leaf curl virus Transcriptional activator protein (C2, L2)

Product Code	CSB-BP329707TGU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27262
Product Type	Recombinant Protein
Immunogen Species	Tomato yellow leaf curl virus (strain Israel) (TYLCV)
Purity	>85% (SDS-PAGE)
Sequence	MQPSSPSTSH CSQVSIKVVH KIAKKKPIRR KRVDLDCGCS YYLHLNCNNH GFTHRGTHHC SSGREWRFYL GDKQSPLFQD NRTQPEAISON EPRHHFHSK IQPQHQEGNG DSQMFSRLPN LDDITASDWS FLKSI
Source	Baculovirus
Target Names	C2
Protein Names	Recommended name: Transcriptional activator protein Short name= TrAP Alternative name(s): Protein C2 Protein L2
Expression Region	1-135
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.