



Recombinant Agropine synthesis reductase (mas1)

Product Code	CSB-MP329763AEW
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27874
Product Type	Recombinant Protein
Immunogen Species	Agrobacterium rhizogenes
Purity	>85% (SDS-PAGE)
Sequence	MPLPRAYLGT AQPPASVKEF YFVRHGATDL NEKEMHLQGE KHWGVQGAGT NIGLNGTGKR QAVLAGNVLR KLPIGSVVCSS PLLRAIQTAL IANIGFLCFD IDEDLKERDF GKHEGGYGPL KMFEDNYPDC EDTEMFSLRV AKALTHAKNE NTLFVSHGGV LRVIAALLGV DLTKEHTNNG RVLHFRRGFS HWTVEIHQSP VILVSGSNRG VGKAI AEDLI AHGYRLSLGA RKVKDLEVAF GPQDEWLHYA RFDAEDHGTM AAWVTA AVEK FGRIDGLVNN AGYGEPVNL D KHV D YQR FHL QWYINC VAPL RMT ELCLPHL YETGSGRIVN INSMSGQRVL NPLVGYNMTK HALGGLTKTT QHVGWDRRCA AIDICLGFVA TDM SAWTDLI ASKDMIQPED IAKLVREAIE RPNRAYVPRS EVMCIKEATR
Source	Mammalian cell
Target Names	mas1
Protein Names	Recommended name: Agropine synthesis reductase EC= 1.-.-.
Expression Region	1-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.