



Recombinant Serine/threonine-protein kinase-transforming protein Rmil (V-RMIL)

Product Code	CSB-BP329771ATD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27966
Product Type	Recombinant Protein
Immunogen Species	Avian rous-associated virus type 1
Purity	>85% (SDS-PAGE)
Sequence	MEAVIKDLIR DQGVREGGGS TAGLSATPPA SLPGLSTNVK ALQKSPGPQR ERKSSSSSED RNRMKTLGRR DSSDDWEIPD GQITVGQRIG SGSFGTVYKG KWHGDVAVKM LNVTAPTPQQ LQAFKNEVGV LRKTRHVNIL LFMGYSTKPQ LAIVTQWCEG SSLYHHLHII ETKFEMIKLI DIARQTAQGM DYHAKSIIH RDLKSNNIFL HEDLTVKIGD FGLATVKSRW SGSHQFEQLS GSILWMAPEV IRMQDKNPYS FQSDVYAFGI VLYELMTGQL PYSNINNRDQ IIFMVGRGYL SPDLKVRSN CPKAMKRLMA ECLKKKRDER PLFPQILASI ELLARSLPKI HRSASEPSLN RAGFQTEDFS LYACASPKTP IQAGGIGEWA VHLLKGLLLG LVVILLVVC LPCLLQCVSS SIRKIDNSL GYREEYKKLQ EAYKQPERRA
Source	Baculovirus
Target Names	V-RMIL
Protein Names	Recommended name: Serine/threonine-protein kinase-transforming protein Rmil EC= 2.7.11.1
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.