



Recombinant Arabidopsis thaliana TATA-box-binding protein 1 (TBP1)

Product Code	CSB-MP329782DOA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P28147
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTDQGLEGSN PVDLSKHPSG IVPTLQNIIVS TVNLDCKLDL KAIALQARNA EYNPKRFAAV IMRIREPKTT ALIFASGKMV CTGAKSEDFS KMAARKYARI VQKLGFPKAF KDFKIQNIVG SCDVKFPIRL EGLAYSHAAF SSYEPELFPG LIYRMKVPKI VLLIFVSGKI VITGAKMRDE TYKAFENIYP VLSEFRKIQQ
Source	Mammalian cell
Target Names	TBP1
Protein Names	Recommended name: TATA-box-binding protein 1 Alternative name(s): TATA sequence-binding protein 1 Short name= TBP-1 TATA-binding factor 1 TATA-box factor 1 Transcription initiation factor TFIID TBP-1 subunit
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.