



Recombinant Simian foamy virus type 1 Protein Bel-1 (bel1)

Product Code	CSB-MP329880SGP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29169
Product Type	Recombinant Protein
Immunogen Species	Simian foamy virus type 1 (SFVmac) (SFV-1)
Purity	>85% (SDS-PAGE)
Sequence	MASWEAQEEL RELLHHPED DPPADLTHLL ELDEMEPKVL CGENPGDEKL KKQVIKTPPM HPSTVTWHFG YKQKEDQQDN IKMRDWVPDP SKMSKSTCKR LILLGLYQAC KAQEIIKMDY DVHWEKSVVN EQYFEVEYNC KMCRTVLHEP MPIMYDPETE LWVKPGRLRG PLGSAVYTLK KHYERCLLTL PSLKGTRLPK RRCNPSRRYE TFREHPPTRK RRSKEGIPTD QQPSTSNGDP MALLSGPCGP HSIQPPSCLL QELPKPEVGS PEMAVAMSGG PFWEEVYGDS IFATPLGSSE DQLLSQFD
Source	Mammalian cell
Target Names	bel1
Protein Names	Recommended name: Protein Bel-1 Alternative name(s): Transactivator of spumavirus Short name= Tas Transcriptional transactivator
Expression Region	1-308
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.