



Recombinant *Drosophila virilis* E3 ubiquitin-protein ligase sina (sina)

Product Code	CSB-EP329898DMP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29304
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila virilis</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSNKINPKRR EPTVAAAVAA ATAVVATNTS SSTGSSAGNT SSANTSSSSS SSLSSAGGGD AGMSADLTSL FECVCFDYV LPPILQCSSG HLVCVSCRSK LTCCPTCRGP LANIRNLAME EVASNVKFPK KHSGYGCTAS LUYTEKTEHE ETCECRPYLC PCPGASCKWQ GPLDLVMQHL MMSHKSITTL QGEDIVFLAT DINLPGAVDW VMMQSCFGHH FMLVLEKQEK YDGHQQFFAI VQLIGSRKEA ENFVYRLELN GNRRLTWEA MPRSIHEGVA SAIHNSDCLV FDTSLAQLFA DNGNLGINVT ISLV
Source	<i>E. coli</i>
Target Names	sina
Protein Names	Recommended name: E3 ubiquitin-protein ligase sina EC= 6.3.2.- Alternative name(s): Seven in absentia protein
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.