



# Recombinant Human Putative fidgetin-like protein 2 (FIGNL2)

<b>Product Code</b>	CSB-MP008678HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	A6NMB9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MHWTPEHAQP LNQWPEQHLD VSSTTPSPAH KLELPPGGRQ RCHYAWAHDD ISALTASNLL KRYAEKYSKV LDSPYERPAL GGYSNASFLN GAKGDPEPWP GPEPPYPLAS LHEGLPGTKS GGGGGSGALG GSPVLGNLNP EPLYAGNACG GPSAAPEYAA GYGGGYLAPG YCAQTGAALP PPPPAALLQP PPPPGYGPSA PLYNYPAGGY AAQPGYGALP PPPGPPAPY LTPGLPAPTP LPAPAPPTAY GFPTAAPGAE SGLSLKRKAA DEGPEGRYRK YAYEPAKAPV ADGASYPAAD NGEGRNGFR AKPPGAAEEA SGKYGGGVPL KVLGSPVYGP QLEPFKFEFPE RAPAPRGGFA VPSGETPKG DPGALELVTS KMVDCGPPVQ WADVAGQGAL KAALEELVW PLLRPPAYPG SLRPPRTVLL FGPRGAGKAL LGRCLATQLG ATLLRLRGAT LAAPGAAEGA RLLQAAFAAA RCRPPSVLLI SELEALLPAR DDGAAAGGAL QVPLLACLDG GCGAGADGVL VVGTTSRPAA LDEATRRRFS LRFYVALPDS PARGQILQRA LAQQGCALSE RELAALVQGT QGFSGGELGQ LCQQAAGAG LPGLQRPLSY KDLEAALAKV GPRASAKELD SFVEWDMYK SGH
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FIGNL2
<b>Protein Names</b>	Recommended name: Putative fidgetin-like protein 2
<b>Expression Region</b>	1-653
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.