



Recombinant Chromosomal replication initiator protein DnaA (dnaA)

Product Code	CSB-BP329910SYN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29440
Product Type	Recombinant Protein
Immunogen Species	Serratia marcescens
Purity	>85% (SDS-PAGE)
Sequence	MSLSLWQQCL ARLQDELPAT EFSMWIRPLQ AELSDNTLAL YAPNRFVLDW VRDKYLNNIN GLLNDFCGTD APLLRFEVGS KPITQVISQT VTASVSSAPA APAARTAAPS RPSWDNAAAQ PELSyrSNVN PKHTFDNFVE GKSNQLARAA ARQVADNPGG AYNPLFLYGG TGLGKTHLLH AVGNIGIMARK ANAKVVYMHS ERFVQDMVKA LQNNAIIEFK RYYRSVDALL IDDIQFFANK ERSQEEFFHT FNALLEGNQQ IILTSDRYPK EINGVEDRLK SRFGWGLTVA IEPPELETRV AILMKKADEN DIRLPGEVAF FIAKRLRSNV RELEGALNRV IANANFTGRA ITIDFVREAL RDLLALQEKL VTIDNIQKTV AEYYKIKVAD LLSKRRRSRV ARPRQMAMAL AKELTNHSLP EIGDAFGGRD HTTVLHACRK IEQLREESH IKEDFSNLIR TLSS
Source	Baculovirus
Target Names	dnaA
Protein Names	Recommended name: Chromosomal replication initiator protein DnaA
Expression Region	1-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.