



Recombinant Escherichia coli Putative uncharacterized protein ychG (ychG)

Product Code	CSB-MP329974ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P30192
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MPLLNDLLDF SDHPLMPPPS AQLFAEHLPT EWIQHCLTLS AHATVRRRRL PGDMVIWMVV QNEPITDVVR RLNLSADGEA GMNLLARSAV TQARQRVGAA PVEWLFQRQA QDRGAERYLK DDWHGLQLFA IDGAQFRTPD KPELREYYGS ANTSTKRQNA YPVMRLVALM NLGSHILLNA VTAPYRQSET VLAHSM LATI PDNSITLFDK LFYSEDLLLT LNQKGCNRHW LLPAWKNIAS EMIELGNTAS PGTIPKRLEH LRGALEVVFI TKRPRPSRPR SVKISKTRYP VKHSAAPLK
Source	Mammalian cell
Target Names	ychG
Protein Names	Recommended name: Putative uncharacterized protein ychG
Expression Region	1-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.