



Recombinant Escherichia coli Protein smf (smf)

| | |
|--------------------------|---|
| Product Code | CSB-YP330032ENV |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P30852 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain K12) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MVDTEIWLRL MSISSLYGDD MVRIAHWVAK QSHIDAVVLQ QTGLTLRQAQ RFLSFPKRSI ESSLCWLEQP NHHLIPADSE FYPPQLLATT DYPGALFVEG ELHALHSFQL AVVGSRAHSW YGERWGRFLFC ETLATRGVTI TSGLARGIDG VAHKAALQVN GVSIAVLGNG LNTIHPRRHA RLAASLLEQG GALVSEFPLD VPPLAYNFPR RNRIISGLSK GVLVVEAALR SGSLVTARCA LEQGREFVAL PGPIGNPGSE GPHWLIKQGA ILVTEPEEIL ENLQFGLHWL PDAPENSFYS PDQQDVALPF PELLANVGDE VTPVDVVAER AGQPVPEVVT QLLELELAGW IAAVPGGYVR LRRACHVVRT NVFV |
| Source | Yeast |
| Target Names | smf |
| Protein Names | Recommended name: Protein smf |
| Expression Region | 1-374 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |