



# Recombinant Escherichia coli 2-dehydro-3-deoxygalactonokinase (dgoK)

<b>Product Code</b>	CSB-YP330084ENV
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P31459
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTARYIAIDW GSTNLRAWLY QGDHCLESRQ SEAGVTRLNG KSPA AVLAEV TTDWREKTP VVMAGMVGSN VGWKVAPYLS VPACFSSIGE QLTSVGDNIW IIPGLCVSHD DNHNVMRGE E TQLIGARALA PSSLYVMPGT HCKWVQADSQ QINDFRTVMT GELHHLLLNH SLIGAGLPPQ ENSADAFTAG LERGLNTPAI LPQLFEVRAS HVLGTLPREQ VSEFLSGLLI GAEVASM RDY VAHQHAILTV AGTSLTARYQ QAFQAMGCDV TAVAGDTAFQ AGIRSIHAV AN
<b>Source</b>	Yeast
<b>Target Names</b>	dgoK
<b>Protein Names</b>	Recommended name: 2-dehydro-3-deoxygalactonokinase EC= 2.7.1.58 Alternative name(s): 2-keto-3-deoxy-galactonokinase 2-oxo-3-deoxygalactonate kinase
<b>Expression Region</b>	1-292
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.