



Recombinant Escherichia coli Putative uncharacterized protein yjbl (yjbl)

Product Code	CSB-MP330189ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32690
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MKKIECACNF LMDKDAQGYI DLSDLTLTSC HFKGDVSKV SFLSSNLQHV TFECKEIGDC NFTAIVDNV IFRCRRLHNV IFIKASGECV DFSKNILDTV DFSQSQLGHS NFRECQIRNS NFDNCYLYAS HFTRAEFLSA KEISFIKSNL TAVMFDYVRM STGNFKDCIT EQLELTIDYS DIFWNEDLDG YINNIKMID TLPDNAMILK SVLAVKLVQM LKILNIVNKN FIENMKKIFS HCPYIKDPII RSYIHSDEDN KFDDFMQRHR FSEVNFDTQQ MIDFINRFNT NKWLIDKNNN FFIQLIDQAL RSTDDMIKAN VWHLYKEWIR SDDVSPIFIE TEDNLRFTNT NELTRNDNIF ILFSSVDDGP VMVVSSQRLH DMLNPTKDTN WNSTYIYKSR HEMLPVNLTQ ETLFSSKSHG KYALFPIFTA SWRAHRIMNK GV
Source	Mammalian cell
Target Names	yjbl
Protein Names	Recommended name: Putative uncharacterized protein yjbl
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.