



Recombinant *Saccharomyces cerevisiae* Serine/threonine-protein phosphatase PP2A-like PPG1 (PPG1)

Product Code	CSB-EP330202SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32838
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MELDECLERL YKAQLLPEVT VRALCFKLKE MLVKESNVIH IQTPVTVVGD MHGQFHDMLE IFQIGGPVPD TNYLFLGDYV DRGLYSVETI MLLIVLKLRY PSRIHLLRGN HESRQITQSY GFYTECLNKY GGNSRVWQYL TDIFDYLVLVLC CIIDDEIFCV HGGLSPNVQT IDQIKIIDRF REIPHDGAMA DLVWSDPEEN NNPTLDHPDN SGQHFQVSPR GAGYTFGRSV VEKFLRMNDM NRIYRAHQLC NEGYQIYFDG LVTTVWSAPN YCYRCGNKAS ILELYSKDQF YFNVFEEAPE NKLLKENS MN DNALEDSISN PVANRKL IAD YFEDDSASAD GSTDPEMYIF SDVYQARSAS NRHVDYFL
Source	<i>E.coli</i>
Target Names	PPG1
Protein Names	Recommended name: Serine/threonine-protein phosphatase PP2A-like PPG1 EC= 3.1.3.16
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.