



Recombinant *Saccharomyces cerevisiae* Protein LTV1 (LTV1)

Product Code	CSB-EP330314SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34078
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	SKKFSSKNS QRYVVVHRPH DDPSFYDTDA SAHVLVPVSN PNKTSPEADL RKKDVSSTKP KGRRAHVGEA ALYGINFDDS EYDYTQHLKP IGLDPENSIF IASKGNEQKV EKKNIEDLFI EPKYRRDEIE KDDALPVFQR GMAKPEYLLH QQDTTDEIRG FKPDMNPALR EVLEALEDEA YVNDVVVE DISKKTQLQG DNYGEEEEKED DIFAQLLGSG EAKDEDEFED EFDEWDIDNV ENFEDENYVK EMAQFDNIEN LEDLENIDYQ ADVRRFQKDN SILEKHNSDD EFSNAGLDSV NPSEEEEDVLG ELPSIQDKSK TGKKKRKSRQ KKGAMSDVSG FSMSSAIAR TETMTVLDDQ YDQIINGYEN YEEEELEED E QNYQPFDSA ERSDFESMLD DFLDNYELES GGRKLAKKDK EIERLKEAAD EVSKGKLSQR RNRERQEKKK LEKVTNTLSS LKF
Source	E.coli
Target Names	LTV1
Protein Names	Recommended name: Protein LTV1 Alternative name(s): Low-temperature viability protein 1
Expression Region	2-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.