



Recombinant Uncharacterized protein D2007.2 (D2007.2)

Product Code	CSB-BP330344CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34376
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MASSSSAALR PFGTARLTPG RQTGRQTQQQ ISAPEPMIAI ALSDMEYLMR SARSSRSSGR SERQPGCTCC CHIAAGVTVI PRVARFSIEG GLSTHTLMNH SDNRIAVKIT CSDNNMYRVT PVYATVEPGQ SLPLHIARIT SDLIKRDRLC VNILEADGNK EAREIFKKNA NTRAPASINM ALEATNDNQN HHHQE
Source	Baculovirus
Target Names	D2007.2
Protein Names	Recommended name: Uncharacterized protein D2007.2
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.