



Recombinant Human Ras association domain-containing protein 10 (RASSF10)

Product Code	CSB-MP019373HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NK89
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDPSEKKISV WICQEEKLVS GLSRRTTCSD VVRVLLEDGC RRRRRQRRSR RLGSAGDPHG PGELPEPPNE DDEDDDEALP QGMLCGPPQC YCIVEKWRGF ERILPNKTRI LRLWAAWGEE QENVRFVLVR SEASLPNAGP RSAEARVLS RERPCPARGA PARPSLAMTQ EKQRRVVRKA FRKLAKLNRR RQQQTPSSCS STSSSTASSC SSSPRTHESA SVERMETLVH LVLSQDHTIR QQVQRLHELD REIDHYEAKV HLDRMRRHGV NYVQDTYLVG AGIELDGSRP GEEPEEVAE AEEAAAAPPL AGEAQAAALE ELARRCDDL RLQEQRVQQE ELLERLSAEI QEELNQRWMR RRQEELAARE EPLEPDGGPD GELLLEQERV RTQLSTSLYI GLRLNTDLEA VKSDLDYSQQ QWDSKKRELQ GLLQTLHTLE LTVAPDGAPG SGSPSREPGP QACADMWVDQ ARGLAKSGPG NDESDTGSL SMHSQDSDSL PMCESLV
Source	Mammalian cell
Target Names	RASSF10
Protein Names	Recommended name: Ras association domain-containing protein 10
Expression Region	1-507
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.