



Recombinant Uncharacterized protein F54H12.5 (F54H12.5), partial

Product Code	CSB-EP330356CXY-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34459
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MLFISRLLP ILYFQLFKKM SPFRLLSLPT LALKNVLLHI DFIDLLELSL ASKKCEIYMK TCCLKIDSLH FHFRRQLNS KSKSLEFDNF SVVDISACSR KVKGSRYDAW TKNYDESR II CIPTNFDEDI IAVFQHFWDL FQIKNLHYDV VGFTGVFLNN TMARLPSK LK SLEFGDYLSG NQDNIDTILN YYDISDNLKV FDEVSKIHEK MKKVNNLSLI CLPSIPIEEL IQLDCETIKL YWKENSAINA VRSLISNWKG NKTRLEHLEI LASEDWNFDL VLEGMNAKPW NPRRRPSNYQ KEDLNIDCTS YMDIEGDCQI ASIGIRENRY LDFIVWF
Source	E.coli
Target Names	F54H12.5
Protein Names	Recommended name: Uncharacterized protein F54H12.5
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.