



# Recombinant Uncharacterized protein F54H12.5 (F54H12.5), partial

<b>Product Code</b>	CSB-MP330356CXY
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P34459
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLFISRSLLP ILYFQLFKKM SPFRLLSLPT LALKNVLLHI DFIDLLELSL ASKKCEIYMK TCCLKIDSLH FHFRRQLNS KSKSLEFDNF SVVDISACSR KVKGSRYDAW TKNYDESR II CIPTNFDEDI IAVFQHFWDL FQIKNLHYDV VGFTGVFLNN TMARLPSK LK SLEFGDYLSG NQDNIDTILN YYDISDNLKV FDEVSKIHEK MKKVNNLSLI CLPSIPIEEL IQLDCETIKL YWKENSAINA VRSLISNWKG NKTRLEHLEI LASEDWNFDL VLEGMNAKPW NPRRRPSNYQ KEDLNIDCTS YMDIEGDCQI ASIGIRENRY LDFIVWF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	F54H12.5
<b>Protein Names</b>	Recommended name: Uncharacterized protein F54H12.5
<b>Expression Region</b>	1-337
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.