



Recombinant Protein kup-1 (kup-1)

Product Code	CSB-EP330387CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P34705
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MDDELDYGSD HSIHGDVRED LSDSELQLPK DRLDAIIDRK EKIMATEKAFG INDSSDKHLN VSIAETAALR KEIGASDQNS HLDSIYAHGV EAMDEFEIQK MFANFRPEKV WKKDNVAMVQ FQYRRDVAAM MLNMSKMMRR VRGRKKADEE GEVLSDDDDV EEGQIMQEKD DDVELIEDLK PNEKGIVASE KNEFITVDI NAREVPNGKW RVLTKHV PAN MFVIIRYATT DEYQTMVTSD RSQVKKTGVK RGNSFWTHE SSNRGGLNVF DKEGKELEWD YEHDFRYNE DKNEEKVEKV KLPQGIKIG RGAVKCGFLF GEGSSSLASD DTTPVKKRRT DEKDYEKDDI VSRQGSSAHA IRPGRVERPI RERIQFPGRE DPDEY
Source	E.coli
Target Names	kup-1
Protein Names	Recommended name: Protein kup-1
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.