



Recombinant *Saccharomyces cerevisiae* Transcription initiation factor TFIID subunit 14 (TAF14)

Product Code	CSB-YP330435SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35189
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVATVKRTIR IKTQQHILPE VPPVENFPVR QWSIEIVLLD DEGKEIPATI FDKVIYHLHP TFANPNRTFT DPPFRIEEQG WGGFPLDISV FLEKAGERK IPHDNLNQLQE SYEVEHVIQI PLNKPLLTEE LAKSGSTEET TANTGTIGKR RTTTNTTAEP KAKRAKTGSA STVKGSDVLE KLAFLGLTKLN EDDLVGVVQM VTDNKTPEMN VTNNVEEGEF IIDLYSLPEG LLKSLWDYVK KNTE
Source	Yeast
Target Names	TAF14
Protein Names	Recommended name: Transcription initiation factor TFIID subunit 14 Alternative name(s): Actin non-complementing mutant 1 Chromosome stability protein 10 SWI/SNF chromatin-remodeling complex subunit TAF14 SWI/SNF complex 29 kDa subunit
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.