



Recombinant Hypoderma lineatum Hypodermin-B

Product Code	CSB-MP330461HRQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35588
Product Type	Recombinant Protein
Immunogen Species	Hypoderma lineatum (Early cattle grub) (Common cattle grub)
Purity	>85% (SDS-PAGE)
Sequence	IVGGFEADIE DFPWQVSIQR GGYHFCGGSI YSPEIIVTAA HCLEKIDASQ LRVRVGGSYW DEEGSLLTVS NFKIHEKYDA MIMWYDVALL KLSSKLTGYA TVKNIELAKE TPPDNADAVV SGWGTIYENY PYMPVQLRSV DVKIVSREVC RSDEYGYGNA IGPTMICAYA VGKDACQGDS GGPLVVGEAL VGVVSWGEGC AYPGFPGVYT DVSVVRSWIT ENAKSF
Source	Mammalian cell
Protein Names	Recommended name: Hypodermin-B Short name= HB EC= 3.4.21.-
Expression Region	31-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.