



Recombinant Cuticle collagen dpy-10 (dpy-10)

Product Code	CSB-MP330489CXY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35800
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	LPIMY NHVQNTITYV DREMAYCEVR FLYFYCIKFF STLQRSNDEA ALELQYGKMR MTGNRTARGA YGSGASHGFR PTAYGDEITG APLETECPGC CIPGPPGPRG SSGTPGKPL PGNAGKPGMP GTTPNQTCPL NQVREPPPCR PCPKGPPGIK GWPGFPGDVG PPGPPGLKGI DGEDGAPGET GPTGPPGYRG GPGAPGDKGP TPEGDLKEGP PGDEGPPGPI GAPGMPGLPG RNGLTGGQGE RGWPGVSGES GEPGYGPEG PMGGQGPPGE PGVCVCQNVD SILLINPGPQ PRIRADDYNS DDGYGGSRRG GDRAGYQGYG RK
Source	Mammalian cell
Target Names	dpy-10
Protein Names	Recommended name: Cuticle collagen dpy-10 Alternative name(s): Protein dumpy-10
Expression Region	46-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.