



Recombinant *Saccharomyces cerevisiae* rRNA-processing protein EBP2 (EBP2)

Product Code	CSB-EP330525SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36049
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAKGFKLKEL LSHQKEIEKA EKLENDLKKK KSQELKKEEP TIVTASNLLK LEKKEKKADV KKEVAADTEE YQSQUALSKKE KRKLKELKK MQEQDATEAQ KHMSGDEDES GDDREEEEEEE EEEEEGRDL EKLAKSDSES EDDSESENDS EEDEDVVAKE ESEEKEEQEE EQDVPLSDVE FSDADVVPK HKLTVNNTKA MKHALERVQL PWKKHSFQEH QSVTSETNTD EHIKDIYDDT ERELAFYKQS LDAVLVARDE LKRLKVPFKR PLDYFAEMVK SDEHMDKIKG KLIEEASDKK AREEARQRQ LKKFGKQVQN ATLQKRQLEK RETLEKIKSL KNKRKHNEID HSEFNVGVEE EVEGKRFDRG RPNGKRAAKN AKYGQGGMKR FKRKNDATSS ADVSGFSSRK MKGKTNRPGK SRRARRF
Source	E.coli
Target Names	EBP2
Protein Names	Recommended name: rRNA-processing protein EBP2 Alternative name(s): EBNA1-binding protein homolog
Expression Region	1-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.