



# Recombinant *Saccharomyces cerevisiae* rRNA-processing protein EBP2 (EBP2)

<b>Product Code</b>	CSB-EP330525SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P36049
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAKGFKLKEL LSHQKEIEKA EKLENDLKKK KSQELKKEEP TIVTASNLKK LEKKEKKADV KKEVAADTEE YQSQUALSKKE KRKLKKELKK MQEQDATEAQ KHMSGDEDES GDDREEEEEEE EEEEEGRDL EKLAKSDSES EDDSESENDS EEDEDVVAKE ESEEKEEQEE EQDVPLSDVE FDSADVVPH HKLTVNNTKA MKHALERVQL PWKKHSFQEH QSVTSETNTD EHIKDIYDDT ERELAFYKQS LDAVLVARDE LKRLKVPFKR PLDYFAEMVK SDEHMDKIKG KLIEEASDKK AREEARQRQ LKKFGKQVQN ATLQKRQLEK RETLEKIKSL KNKRKHNEID HSEFNVGVEE EVEGKRFDRG RPNGKRAAKN AKYGQGGMKR FKRKNDATSS ADVSGFSSRK MKGKTNRPGK SRRARRF
<b>Source</b>	E.coli
<b>Target Names</b>	EBP2
<b>Protein Names</b>	Recommended name: rRNA-processing protein EBP2 Alternative name(s): EBNA1-binding protein homolog
<b>Expression Region</b>	1-427
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.