



Recombinant Escherichia coli Protein transport protein HofB homolog (hofB)

Product Code	CSB-EP330617ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36645
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MNIPQLTALC LRYHGVLLDA SEEVVHVAVV DAPSHELLDA LHFATTKRIE ITCWTRQQME GHASRTQQTL PVAVQEKHQP KAELLTRTLQ SALEQRASDI HIEPADNAYR IRLRIDGVLH PLPDVSPDAG VALTARLKV L GNLDIAEHRL PQDGQFTVEL AGNAVSFRIA TLPCRGGEKV VLRLQLQVVGQ ALDVNTLGMQ PLQLADFAHA LQQPQGLVLV TGPTGSGKTV TLYSALQKLN TADINICSV DPVEIPIAGL NQTQIHPRAG LTFQGVLRAL LRQDPDVIMI GEIRDGETAE IAKAAQTGH LVLSTLHTNS TCETLVRLLQ MGVARWMLSS ALTLVIAQRL VRKLCPHCRR QQGEPIHIPD NVWPSPLPHW QAPGCVHCYH GFYGR TALFE VLPITPVIRQ LISANTDVES LETHARQAGM RTL FENGCLA VEQGLTTFEE LIRVLGMPHG E
Source	E.coli
Target Names	hofB
Protein Names	Recommended name: Protein transport protein HofB homolog
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.