



Recombinant Escherichia coli Uncharacterized protein yccE (yccE)

| | |
|--------------------------|--|
| Product Code | CSB-BP330618ENV |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P36661 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain K12) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MGSNIHGISC TANNYLKQAW NDIKNEYEKN QTYSITLFEN TLVCFMRLYN ELRRKVNEED TPCLECESLE KEFEEMQNDN DLSLFMRILR TNDTQIYSGV SGGITYTIQY VRDIDIVRVS LPGRASESIT DFKGYYWYNF MEYIENINAC DDVFSEYCFD DENISVQPER INTPGISDLD SDIDLSGISF IQRETNQALG LKYAPVDGDG YCLLRAILVL KQHDYSWALV SYKMQKEVYN EFIKMVDKKT IEALVDTAFY NLREDVKTFL GVDLQSDNQi QGQSSLMSWS FLFFKKQFID SCLNNEKCIL HLPEFIFNDN KNLLALD TDT SDRIKAVKNF LVVLSDSICS LFIVNSNVAS ISLGNESFST DEDLEYGYLM NTGNHYDVYL PPELFAQAYK LNNKEMNAQL DYLNRYAI |
| Source | Baculovirus |
| Target Names | yccE |
| Protein Names | Recommended name: Uncharacterized protein yccE Alternative name(s): ORF-D |
| Expression Region | 1-418 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |