



Recombinant Escherichia coli Ribosomal RNA large subunit methyltransferase A (rlmA)

Product Code	CSB-MP330668ENV
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P36999
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSFSCPLCHQ PLSREKNSYI CPQRHQFDMA KEGYVNULLPV QHKRSRDPGD SAEMMQARRA FLDAGHYQPL RDAIVAQLRE RLDDKATAVL DIGCGEGYYT HAFADALPEI TTFGLDVSKV AIKAAAKRYP QVTFCVASSH RLPFSDTSM AIIRIYAPCK AEELARVVKP GGWVITATPG PRHLMELKGL IYNEVHLHAP HAEQLEGFTL QQSAELCYPM RLRGDEAVAL LQMTPFRAWRA KPEVWQTLAA KEVFDCQTFD NIHLWQRSY
Source	Mammalian cell
Target Names	rlmA
Protein Names	Recommended name: Ribosomal RNA large subunit methyltransferase A EC= 2.1.1.187 Alternative name(s): 23S rRNA m1G745 methyltransferase rRNA (guanine-N(1)-)-methyltransferase
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.