



Recombinant Human Ankyrin repeat and SOCS box protein 14 (ASB14)

Product Code	CSB-MP002183HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6NK59
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDNYTSDEDI DEDFDTQLII QQSLQDIYKP GTAQHAPKDE SLHSFLSADY KKIVETIEKV GKEDALSHLT KYHSAFGEAD EIGWIPLHKA AVQLNRKILE ITLSASDPSL WEQTTHNGET PLFLAVSSCL LENATFLLLN GCNPNNAKNFE GNSPLLA AVL RDCYDMAALL INYGADVNLN CANERTALHE AAKLGREDMV KLMLVSGAHP DPQSTYGFTP LALAAQSGHT EIMEMLLRKG KIFCLASDSS SILLEAASGG NPDAVALLLE YGADANIPKN SGHLPIHVAA DRGHLLALKI LIPVTDLAAI KQSGISPVHC AAAGAHPQCL ELLIQAGFDV NFMLDQRINK HYDDHRKSAL YFAVNSDLS SVKLLLSAGA LPNQDPVNCL QIALRMGNYE LISLLLRHGA NVNYFCRVNP LHFPSALQYT LKDEVMLRML LNYGYDTERC FDCPHGDKVH PSYTVEGWTS TVIKDTKFCE VITLSWLQHL SGKVVRVMLD YVDQVRICKS LKAVLQKQGI WSEIHFILTN PRSLKHLCLRL KIRKCMGRLH LRCPVFMSFL PLPNRLKAYV LYKEYDLYGQ GIFTGTW
Source	Mammalian cell
Target Names	ASB14
Protein Names	Recommended name: Ankyrin repeat and SOCS box protein 14 Short name= ASB-14
Expression Region	1-587
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.