



Recombinant *Saccharomyces cerevisiae* UPF0657 nucleolar protein YBR141C (YBR141C)

Product Code	CSB-MP330852SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38278
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MHSRKSKSIT GKRKQVGSNV TRVIKPQKTR RIIRRFHHLI NKRQSICKFL CLKENLDDSN EEKNDKIIRL SIKGNVRLGK YYEDGKSQSF NDAMESQLLR LHSLIKNESK SKDTSDLAVM YTLLGYIMNQ INKLGGETY QIASQNGQLK ERGGDTSKLL EKWIRSSFEN CPGAVALEIG SLSSGNRISR CALFRNVVRI DLEEHEGVK QDFMERPLPR NENDKFDLIS CSLVLNFVKN HRDRGAMCHR MVKFLKPQGY IFIVLPQACV THSRYCDKTL LQNLLGSIGL IMLNSHQSNK LYYCLYQLQV VPPQPSSFSK RIKVNDGPGL NNFGITL
Source	Mammalian cell
Target Names	BMT2
Protein Names	Recommended name: UPF0657 nucleolar protein YBR141C
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.