



Recombinant Copper-containing nitrite reductase (nirK)

Product Code	CSB-MP330892AGN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38501
Product Type	Recombinant Protein
Immunogen Species	Alcaligenes faecalis
Purity	≥85% (SDS-PAGE)
Sequence	QGA VRKA TAAEIAALPR QKVELVDPPF VHAHSQVAEG GPKVVEFTMV IEEKKIVIDD AGTEVHAMAF NGTVPGPLMV VHQDDYLELT LINPETNTLM HNIDFHAATG ALGGGGLTEI NPGEKTILRF KATKPGVVFVY HCAPPGMVPW HVVSGMNGAI MVLPREGLHD GKGKALTYDK IYYVGEQDFY VPRDENGKYK KYEAPGDAYE DTVKVMRTL TPTHVVFNGAV GALTGDKAMT AAVGEKVLIV HSQANRDTRP HLIGGHGDYV WATGKFNTTP DVDQETWFIP GGAAGAAFYT FQQPGIYAYV NHNLIEAFEL GAAAHFKVTG EWNDDLMTSV LAPSGT
Source	Mammalian cell
Target Names	nirK
Protein Names	Recommended name: Copper-containing nitrite reductase EC= 1.7.2.1 Alternative name(s): Cu-NIR
Expression Region	34-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.