



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YHR054C (YHR054C)

Product Code	CSB-EP330929SVG-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38780
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVPAENLSP IPASIDTNDI PLIANDLKLL ETQAKLINIL QGVPFYLPVN LTKIESLLET LTMGVSNTVD LYFHDNEVRK EWKDTLNFN TIVYTNFFLF VQNESSLSMA VQHSSNNKT SNSERCAKDL MKIISNMHIF YSITFNFIFF IKSIKSFSSG NNRFHSGKE FLFANHFIEI LQNFIAITFA IFQRCEVILY DEFYKNLSNE EINVQLLLIH DKILEILKKI EIIVSFLRDE MNSNGSFKSI KGFNKVLNLI KYMLRFSKKK QNFARNSDNN NVDYSQSAK NKNVLLKFPV SELNRIYLKF KEISDFLMER EVVQRSIID KDLESDNLGI TTANFNDFYD AFYN
Source	E.coli
Target Names	YHR054C
Protein Names	Recommended name: Uncharacterized protein YHR054C
Expression Region	1-354
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.