



Recombinant *Saccharomyces cerevisiae* Ribosome biogenesis protein SSF1 (SSF1)

Product Code	CSB-BP330931SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38789
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MAKRRQKKRT HAQLTPEQEQ GIPKSMVIRV GQTS LANHSL NQLVKDFRQI MQPHTAIK LK ERKSNK LKDF VVMCGPLGVT HLFMFTQSEK TGNVSLKIAR TPQGPTVTFQ VLDYSLGRDI KKFLKRPKSL NNDDVLNPPL LVLNGFSTSK RSGEDDQDVN VEKVIVSMFQ NIFPPLNPAR TSLNSIKRVF MINKDRETGE ISMRHYFIDI REVEISRNLK RLYKAKNNLS KTVPNLHRKE DISSLILDHD LGAYTSESEI EDDAIVRVVD NQDVKAKHSQ SLKSQRTPVE KKD NKEREKE TEEEDVEMEE PKPSEN LQPT PRKKA IKLTE LGPRLTLKLV KIEEGICSGK VLHHEFVQKS SEEIKALEKR HAAKMRLKEQ RKKEQEENIA KKKAVKDAKK QRKLERRKAR AAEGGEGQ GK DDAMS DDESS SSDSEHYGSV PEDLDSDLFS EVE
Source	Baculovirus
Target Names	SSF1
Protein Names	Recommended name: Ribosome biogenesis protein SSF1
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.