



# Recombinant *Saccharomyces cerevisiae* tRNA (guanine (37)-N1)-methyltransferase

<b>Product Code</b>	CSB-BP330932SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P38793
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	KIPMCA VKFPEPKNIS VFSKNFKNCI LRVPRIPHVV KLNSSKPKDE LTSVQNKKLK TADGNNTPTV KGVLLHESIH SVEDAYGKLP EDALAFLEN SAEIVPHEYV LDYDFWKAEE ILRAVLPEQF LEEVPTGFTI TGHIAHLNLR TEFKPFDSL I GQVILDKNNK IECVVDKVSS IATQFRFTFPM KVIAGKSDSL VVEQKESNCT FKDFDSKVYV NSRLHTEHER LVKQYFQPGQ VVCDVFAGVG PFAVPAGKKD VIVLANDLNP ESYKYLKENI ALNKVAKTVK SFNMDGADFI RQSPQLLQQW IQDEEGGKIT IPLPLKKRHR SQQHNDQQPP QPRTKELIIP SHISHYVMNL PDSAISFLGN FRGIFAAHTK GATDTIQMPW VHVHCFEKYP PGDQVTEDEL HARVHARIA ALKVTADDLP LNAVSLHLVR KVAPTKPMYC ASFQLPANV
<b>Source</b>	Baculovirus
<b>Target Names</b>	TRM5
<b>Protein Names</b>	Recommended name: tRNA (guanine(37)-N1)-methyltransferase EC= 2.1.1.228 Alternative name(s): M1G-methyltransferase tRNA [GM37] methyltransferase tRNA methyltransferase 5
<b>Expression Region</b>	45-499
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.