



Recombinant *Saccharomyces cerevisiae* E3 ubiquitin-protein ligase DMA1 (DMA1)

Product Code	CSB-MP330936SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38823
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSTNTVPSSP PNQTPPAASG IATSHDHTKF NNPIRLPISI SLTINDTPNN NSNNNSVSNG LGILPSRTAT SLVVANNGSA NGNVGATAAA AATVETNTAP AVNTTKSIRH FIYPPNQVNQ TEFSLDIHLP PNTSLPERID QSTLKRRMDK HGLFSIRLTP FIDTSSTSVA NQGLFFDPII RTAGAGSQII IGRYTERVRE AISKIPDQYH PVVFKSKVIS RTHGCFKVDD QGNWFLKDVK SSSGTFLNHQ RLSSASTTSK DYLLHDGDII QLGMDFRGGT EEIYRCVKMK IELNKSWKLLK ANAFNKEALS RIKNLQKLTT GLEQEDCSIC LNKIKPCQAI FISPCAHSWH FHCVRRLVIM NYPQFMCPNC RTNCDLETTL ESESESEFEN EDEDEPDIEM DIDMEINNNL GVRLVD
Source	Mammalian cell
Target Names	DMA1
Protein Names	Recommended name: E3 ubiquitin-protein ligase DMA1 EC= 6.3.2.- Alternative name(s): Checkpoint forkhead associated with RING domains-containing protein 1 Defective in mitotic arrest protein 1
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.