



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YHR127W (YHR127W)

Product Code	CSB-YP330938SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38833
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MARNRTTSKK NVQSKRLIDR VVPMDKIKKV GVAKKKTVEH TKEGFSVVNG KLVSSNDVGV LLREAQGAID KRTNVSQRNR KKGIKNNRPH KDINSSPDWG NAHRGTDWQS EKANGMNRK NSRNFTTNIK LQRQHFGEEI QGGSQLVIST NSDASDKLLM LFNLTGVTNQ ENLKNVLENI SQVQIAQIRV RDLPSGSATA KVRLAYPTTQ SLEKVRKLFH GALVDGRRRIQ VVIASDESSH LSY
Source	Yeast
Target Names	YHR127W
Protein Names	Recommended name: Uncharacterized protein YHR127W
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.