



Recombinant *Saccharomyces cerevisiae* Saccharopine dehydrogenase [NAD (+), L-lysine- forming]

Product Code	CSB-BP330963SVG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P38998
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	AAVTLHLRA ETKPLEARAA LPTTVKCLI AKGFKIYVED SPQSTFNINE YRQAGAIIVP AGSWKTAPRD RIIIGLKEMP ETDTFPLVHE HIQFAHCYKD QAGWQNVLMR FIKGHGTLYD LEFLENDQGR RVAAFGFYAG FAGAALGVRD WAFKQTHSDD EDLPAVSPYP NEKALVKDVT KDYKEALATG ARKPTVLIIG ALGRCSGSAI DLLHKVGIPD ANILKWDIKE TSRGGPFDEI PQADIFINCI YLSKPIAPFT NMEKLNNPNR RLRTVVDVSA DTTNPHNPIP IYTVATVFNK PTVLVPTTAG PKLSVISIDH LPSLLPREAS EFFSHDLLPS LELLPQRKTA PVWVRAKKLF DRHCARVKRS SRL
Source	Baculovirus
Target Names	LYS1
Protein Names	Recommended name: Saccharopine dehydrogenase [NAD(+), L-lysine-forming] Short name= SDH EC= 1.5.1.7 Alternative name(s): Lysine--2-oxoglutarate reductase
Expression Region	2-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.