



# Recombinant *Saccharomyces cerevisiae* Poly (A) ribonuclease POP2

<b>Product Code</b>	CSB-BP330965SVG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P39008
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MQSMNVQPRV LAVGGEQFFS QRQASEQHQQ QNMGPQVYSP KVNRRARMFPQ GMPVNTINGS VNQEMNNAYL LKQKNEPLLT QQQQQQQQQQ QPFNIGTPVS VASLPPGLNV LQQQQQQQQQ QQQQQQGVGL NRPLASQLPK HLTNQSMPPI FLPPPNYLFV RDVWKSNNLYS EFAVIRQLVS QYNHVSISTE FVGTLARPIG TFRSKVDYHY QTMRANVDFL NPIQLGLSLS DANGNKPDNG PSTWQFNFEF DPKKEIMSTE SLELLRKSGI NFEKHENLGI DVFEFSQLLM DSGLMMDDSV TWITYHAAYD LGFLINILMN DSMPNNKEDF EWWWVHQYMPN FYDLNLVYKI IQEFKNPQLQ QSSQQQQQQQ YSLTTLADEL GLPRFSIFTT TGGQSLLMLL SFCQLSKLSM HKFPNGTDFA KYQGVIIYGID GDQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	POP2
<b>Protein Names</b>	Recommended name: Poly(A) ribonuclease POP2 EC= 3.1.13.4 Alternative name(s): CCR4-associated factor 1
<b>Expression Region</b>	1-433
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.